




SOUTHERNIMPLANTS®


Innovative Treatment Solutions

SOFT BONE Implants
Product Catalogue



Ø4.0mm Implants (Tapered)

Tapered,
Soft Bone



PT-DC40xx

(where **xx** is implant length)

Implants are pre-mounted and available in lengths of:

NOTE: Implant dimensions and information - page 13.

| ITEM CODE | IMPLANT LENGTHS (in mm) |
|-----------|-------------------------|
| PT-DC40 | 13 / 15 / 18 / 20 |



Surgical Components

Cover Screw

CS-DC4



Healing Abutments

HA-DC4-N2

Ø4.0



HA-DC4

Ø4.5



4/6
lengths

HA-DC4-W

Ø5.5



2/4/6
lengths

Prosthetic Flowchart

DIRECT

Healing Abutments

HA-DC4-N2



HA-DC4



HA-DC4-W



Impression Copings

IC-DC4
(pick-up)



IC-DC4-W
(pick-up)



SFT-DC4
(scanning flag)



Titanium

ICT-DC4
(transfer)



ICT-DC4-W
(transfer)



Laboratory Analogues

LA-DC4



LAD-DC4



(digital analogue)

Prosthetic Components

GC-DC4
(engaging)
GC-NDC4
(non-engaging)



Gold

TC-DC4-1
(engaging)
TC-NDC4-1
(non-engaging)



Titanium

PKC-DC4-2
(engaging)
PKC-NDC4-2
(non-engaging)



PEEK*

*Screw Torque: 20Ncm

TCA-DC4
(engaging)



TCA12-DC4
(engaging)



TCA24-DC4
(engaging)



CAB-DC4
(engaging)



TIB-DC4
TIB-DC4-C1.5
TIB-DC4-C3
(engaging)



TIB-NDC4
TIB-NDC4-C1.5
TIB-NDC4-C3
(non-engaging)



TIB abutments available in 3 cuff height variants:
0.7mm, 1.5mm and 3mm

PA-DC4
(engaging)
PA-NDC4
(non-engaging)



Passive Abutment Screws

PA-DC4-16T



Titanium

TORQUE:
15Ncm

PA-DC4-16B



(Anodized blue
and for laboratory
use only)

TS-DC4-16



Titanium

TORQUE:
20Ncm

BS-DC4-16



Brass*

*(Blackened and
for laboratory
use only)

INDIRECT

Compact Conical Abutments

MC-DC4



1 / 3 / 4 / 5

MC-DC4-20D*



MC-DC4-30D*



*Packed with
TS-DC4-MC

HMC



4/6

OR



4/6

CMC1
(pick-up)



CMC2
(transfer)



CMC-ZG-2
(transfer)



LSMC1



LAD-MC



(digital analogue)

GMC1



Gold

TMC1 / 5



Titanium
TMCSL
(long version)

PKC-MC



PEEK*

*Screw Torque:
10 - 15Ncm

TIB-MC-48
(titanium
scanning abutment)



PA-MC-48
(passive abutment)



1 Series
Screws

TORQUE:
10 - 15Ncm


NOTE: Refer to CAT-1189 for Equator Overdenture Abutments.

Site Preparation Sequence

Ø4.0mm Tapered (PT-DC40xx)

(Illustration is for 13mm implant)


Soft Bone



NOTE: Site preparation sequence recommended by Southern Implants does not replace the judgement and experience of the surgeon.

IMPLANT DIMENSIONS AND INFORMATION

| RANGE | | MAJOR DIAMETER | PLATFORM WIDTH | PROSTHETICS | CONNECTION DIAMETER | THREAD PITCH | APEX DIAMETER | Cylindrical or Tapered | PLATFORM ANGLE | IMPLANT LENGTHS | | | |
|-----------|--------|----------------|----------------|-------------|---------------------|--------------|---------------|------------------------|----------------|-----------------|----|----|----|
| | | | | | | | | | | 13 | 15 | 18 | 20 |
| ● PT-DC40 | Ø4.0mm | 4.0 | 3.70 | ● | 2.80 | 1.6 | 2.0 | T | - | ✓ | ✓ | ✓ | ✓ |



All dimensions in this catalogue are in mm, unless otherwise specified.
Not all products are cleared for sale in all countries.

For more information, please contact your
Southern Implants Representative or visit southernimplants.com



South Africa - Headquarters 1 Albert Road, Irene, RSA
T: +27-12-667-1046 | **E:** info@southernimplants.com

EC REP
Southern Implants Europe AB Holmgatan 30, S-791 71, Falun, Sweden
E: ecrep@southernimplants.com

Subsidiaries

Australia
Southern Implants Australia
T: +61-2-8076-9337
E: info@southernimplants.com.au

Spain and Portugal
Southern Implants Ibérica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland
Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-099
E: customercare@southernimplants.com

Southern Implants are distributed world-wide, please visit southernimplants.com for a list of Distributors.